LIFEWAVE®

# SP6 Complete

## **Instructions**

Mode d'emploi | Instrucciones | Instruksi

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## INSTRUCTIONS FOR USE


Place one SP6 Complete patch on the body, using any of the recommended locations shown on page 8. For optimal results, rotate the patch between the placements that provide you the most control of your appetite.

Apply the patches to clean, dry skin in the morning. Patches may be worn for up to 12 hours before discarding. Keep well hydrated while using this product. Do not use more than one SP6 Complete patch at a time on the body.

### STEP 1

First, locate the acupressure point you wish to stimulate.

### STEP 2

Next, remove a plastic bead from the bag. Remove the backing from the adhesive plastic patch included. Place the plastic bead in the center of the patch on the adhesive side of the patch.

### STEP 3

Lastly, apply the bead and patch to the acupressure point, smoothing the patch on so that it is properly secured to the skin. There should be a mild pressure on the acupressure point.

**Warnings:** Remove immediately if you feel discomfort or skin irritation occurs. Do not reuse patch once removed from the skin. For external use only. Do not ingest. Do not use on wounds or damaged skin. Ask a health professional before using if you have a health condition, any questions or concerns about your health. Do not use if pregnant or nursing.

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